


# A fresh start for dental chairs

DÜRR SYSTEM HYGIENE




Cleaning and care for vinyl upholstery with FD 360

**DÜRR**  
SYSTEM-  
HYGIENE



Even heavily soiled vinyl upholstery can be visibly restored and deep-cleaned with FD 360



Dürr FD 360 – 500 ml in the pump spray bottle in a set with 10 special sponges

## Restoration and maintenance of value for highly stressed vinyl upholstery in the surgery – Dürr FD 360

Dental chairs are the centre of every surgery. They are highly stressed and require regular cleaning and care to maintain their economic value in the long run. The new FD 360 cleaning and care for vinyl upholstery can be applied easily and cleans thoroughly:

- ▶ Intensive cleaning effect – removes dirt and stains of all kind, dirt rims and textile discolouring
- ▶ Outstanding compatibility with all synthetic leathers, tested on treatment units by leading equipment manufacturers
- ▶ Practical set containing 500-ml spray bottle and 10 special sponges for a particularly effective application

**DÜRR**  
DENTAL

# Dürr FD 360 cleaning and care for vinyl upholstery

## DÜRR SURFACE DISINFECTION AND CLEANING

### Perfectly clean and thorough

The dental chair is a highly stressed equipment in a inventory. As part of the furniture it instantly catches the patient's eye, which is why its condition is the hallmark of a surgery. Immaculate hygiene and care are therefore very important. However, in the course of time the state of the upholstery can degrade visibly through regular wear and tear. Waterspots of the instrument cooling water spraymist, patients' sweat, as well as dirt rims and discolouring from textiles leave their ugly traces. With FD 360 unsightly vinyl upholstery can be restored fast and thoroughly for better long-term value.


Because the cleaning agent removes all stains quickly – even intensive, non-fixed colouring agents.

### Maintenance of value and care


While cleaning vinyl upholstery, FD 360 is very gentle on the materials. To remove dirt from stressed surfaces with a long-lasting effect first, the cleaning components make them spotlessly clean; then the care components produce a shiny finish. Thus a heavily soiled dental chair can be completely restored with the first treatment with FD 360 and then maintained in an immaculate condition with only a few millilitres per week.

### FD 360 cleaning for vinyl upholstery

The odourless cleaning fluid FD 360 provides for the deep-cleaning and care for all vinyl upholstery surfaces in the surgery. The emulsion is based on a mix of special tensides, silicon particles and nurturing avocado oil. The fluid is used undiluted and is easily applied. Used with the special application sponges it removes dirt and stains of all kind without a trace – producing a shine. FD 360 is supplied in an 500-ml pump spray bottle with 10 special application sponges.



*In an environmentally friendly procedure, FD 360 is applied from a spray bottle and wiped with a special sponge.*



*Intensive cleaning power with the best material compatibility – this is what characterises FD 360.*



*Vinyl upholstery of all kind and manufacturers benefit from care with FD 360 that leaves a silky shine.*

Dürr Dental GmbH & Co. KG  
Höfigheimer Strasse 17  
D-74321 Bietigheim-Bissingen

[www.duerr.de](http://www.duerr.de)  
[info@duerr.de](mailto:info@duerr.de)

